



# PCRRun®

## Canine Leishmania Molecular Detection Kit

Cat. No.30CLI108

For *in vitro* veterinarian diagnostic use only

User Manual

### INTENDED USE

PCRRun® Canine Leishmania Molecular Detection Kit is intended for detection of *Leishmania infantum* in **DNA** isolated from conjunctival swabs, lymph node fine needle aspirations (FNA), bone marrow and blood using an appropriate DNA extraction kit.

The kit should be employed during an active infection. It contains all the disposable components required for performing an easy and accurate test.

### PRINCIPLE

PCRRun® is a molecular assay based on isothermal amplification of part of kinetoplastid DNA associated with *L. infantum*. It is intended for the qualitative detection of *Leishmania infantum*. This kit is designed to be used with a compatible heat block.

### STORAGE AND HANDLING

- Storage 2-25°C (room temperature or refrigerated).
- Avoid exposure to direct sunlight.
- Do not use beyond expiration date stated on the package label.
- Do not freeze!

#### Precautions:

- The PCRRun® assay is not to be used on the specimen directly. An appropriate nucleic acid extraction method must be conducted prior to using this kit.
- In order to avoid the introduction of contaminating DNA into the reaction, clean laboratory examination gloves must be worn while performing the test.
- Open and remove the PCRRun® reaction tubes from the sealed pouches only immediately prior to their use.
- **Return unused PCRRun® reaction tubes to the original aluminum packet together with the desiccator and seal the zip-lock tightly.**
- Do not use kit if the pouch or the components are damaged.
- Each component in this kit is suitable for use only with a specific lot number. Components have been quality control approved as a standard batch unit. Do not mix components from different lot numbers.
- Employ accepted sanitary procedures designated for biological and molecular waste when handling and disposing of the reaction tubes.

### BACKGROUND

Leishmaniasis is an important vector-borne disease resulting from a protozoan infection caused by parasites classified in the family *Trypanosomatidae*. *Leishmania* spp. are usually transmitted indirectly between hosts by sandflies which act as biological vectors. Mammals can be infected asymptotically for long periods, and they often remain chronically infected even after clinical cure. Subclinically infected animals can transmit *Leishmania* to sandflies and the parasites can also be transmitted via blood transfusions and transplacental transmission.

*L. infantum* is often classified as a visceral disease in humans, however dogs usually have both visceral and cutaneous involvement. Canids have been reported to be the major reservoir hosts for *L. infantum*, and dogs are the most important species in maintaining this parasite in domestic cycles. The reported incubation period for *L. infantum* in

dogs varies from three months to seven years. In some dogs, severe clinical signs occur soon after the animal becomes infected while other dogs remain asymptomatic, in some cases for a lifetime. These animals can demonstrate clinical signs at any time, particularly if they become immunosuppressed. The disease is usually slowly progressive.

### DIAGNOSIS

Leishmaniasis may be manifested as a subclinical infection, self-limiting disease or the classical non-self-limiting severe illness. The main clinical findings associated with typical canine leishmaniasis are skin lesions. Variable symptoms such as lethargy, anorexia, weight loss, anemia, splenomegaly and local or general lymphadenopathy can accompany such an infection. Chronic renal disease, resulting in death, is common with *L. infantum* infections. Erosive or nonerosive polyarthritis with polymyositis can lead to progressive muscle atrophy<sup>3</sup>.

Leishmaniasis can be diagnosed by direct microscopic observation of amastigotes in stained macrophage present in blood, ocular granulomas, tissue collected from lymph node and bone marrow aspirates as well as skin scraping from lesions. Parasites are sometimes undetectable by this method due to low numbers. *Leishmania* spp. can be cultured in a variety of media however *in vitro* culture requires 5-30 days before diagnosis can be determined. Serological testing is available in the form of Indirect Fluorescent Antibody Test (IFAT), Enzyme-Linked Immunosorbent Assay (ELISA) and various agglutination methods. Most symptomatically infected dogs are seropositive, however not all asymptotically infected animals or those with localized skin lesions have detectable antibody titers. Polymerase chain reaction (PCR) has been shown to be more sensitive than conventional techniques of microscopy and culture used alone or in combination<sup>1,2</sup>. PCRRun® is a sensitive method to detect *Leishmania*.

### KIT CONTENTS

Components	Contents	Amount
Aluminum pouch Cat No. 03CLI100	PCRRun® strip of 8 lyophilized <i>Leishmania</i> single reaction tubes	1
PCRRun® buffer Cat No. 31000120	PCRRun® buffer to re-dissolve lyophilized reaction pellets	1 Vial, 200 µl
Detection device Cat No. 03100010	Aluminium pouch with disposable nucleic acid detection device.	8
Capillary tubes Cat No. 03200005	Disposable plastic capillary tubes - 5 µl*	10
Capillary tubes Cat No. 03200015	Disposable plastic capillary tubes - 15 µl*	10

\*Accurate laboratory pipettes with aerosol barrier tips can be used in place of the plastic pipettes.

#### EQUIPMENT TO BE SUPPLIED BY USER:

- PCRRun Sample Prep DNA Extraction Kit is Recommended (Cat No. 30PRE208).
- Heat block which maintains 60°C – compatible with 0.2 PCR tubes. Heat block can be supplied by Biogal
- Timer
- Small scissors
- Fine tipped permanent marker
- Protective laboratory gloves

### SAMPLE COLLECTION, STORAGE AND TRANSPORT

The kit is suitable for detecting nucleic acid extracted from conjunctival swabs, fine needle aspirations (FNA), bone marrow and blood using the above mentioned extraction kit.

For optimal results, recently acquired samples and freshly prepared DNA extracts are recommended for use with the PCRRun® kit.

Unless otherwise recommended, samples and DNA extracts can be maintained at 2-8°C for up to 24 hours or -20°C for an extended period of time.

## PROTOCOL - PCRun® REACTION

1. Prepare a clean working area for the assay. Work area should be decontaminated with commercial household bleach (3.5%) diluted 1:10 with water.
2. Prepare all parts of the assay:

- ✓ Extracted DNA sample
- ✓ Pouch with reaction tubes
- ✓ Capillary tubes for dispensing 5 µl and 15 µl volume
- ✓ Fine tipped permanent marker

3. Switch on the heat block and adjust to 60°C. Once the block has reached the target temperature, continue with the reaction.
4. Remove the PCRun® strip from its protective pouch. Take care to return the unused tubes to the aluminium envelope and seal completely with zip-lock to maintain a dry environment. Eight individual reaction tubes are connected by a thin plastic spacer. Employing a small clean scissors, disconnect the required number of tubes without disturbing the lids. Tap the tubes lightly on a surface and observe that the small white pellet is located on the bottom of the tube.
5. Label the lid of the tubes clearly for sample identification.

6. Carefully open the lid of the reaction tubes, one at a time. Employing the 15 µl disposable capillary tube, dispense 15 µl PCRun® Buffer into the reaction tube. Make sure that the entire content of the capillary tube has been emptied into the PCRun® reaction tube. Tap the tube on a surface to bring all the fluid to the bottom of the tube. Incubate the mixture at room temperature for 1 minute to allow the pellet to completely dissolve. Using the 5 µl capillary tube dispense 5 µl of extracted DNA into the same tube. Mix by lightly tapping the tube. Make sure all fluid is at the bottom of the reaction tube.



7. Place the reaction tube into the appropriate hole in the pre heated block (60°C) and incubate for exactly 1 hour. Do not open the tube cover during or at the end of the incubation period.

8. At the end of the incubation period (1 hr) remove the tube from the heat block and analyze immediately with the disposable nucleic acid detection device.

## ANALYSIS OF PCRun® REACTION WITH THE DISPOSABLE NUCLEIC ACID DETECTION DEVICE

One disposable nucleic acid detection device is needed for each test. Open and remove the components of the detection device. The device consists of two plastic parts, the Amplicon Cartridge containing a plastic buffer bulb and the Detection Chamber containing the lateral flow strip (Figure 1).


1. Verify the presence of fluid in the bulb.
2. Mark each chamber with the sample ID.
3. Align the lid section of the PCRun® reaction tube with the wide partition located beside the buffer bulb. Apply light pressure to attach the reaction tube to the Amplicon Cartridge (Figure 1).
4. Fold the Amplicon Cartridge in two and snap closed. Place the cartridge into the Detection Chamber with the bulb facing downwards and away from the chamber lever.
5. Push the lever downwards to lock the device.
6. Wait for 15-30 minutes to read the results. Results read after 30 minutes are invalid.



## READING AND INTERPRETING THE RESULTS

A valid test must present a red control band. The control line must appear regardless of a positive or negative result. (Figure 2):

1. **Positive Result** - two bands appear, the upper control line and the lower test line. The appearance of both control line and test line indicates the presence of *Leishmania infantum* DNA.
2. **Negative Result** - a single control line appears. The appearance of a control line only, indicates the absence of the *Leishmania infantum* DNA or that the copy number is below the detection limit.



## LIMITATIONS

As with any diagnostic kit, the results obtained must be used as an adjunct to other clinical and laboratory findings. The accuracy of the test results depends on the quality of the sample and adherence to protocol. A negative result may be obtained if the specimen is inadequate or target pathogen concentration is below the sensitivity of the test.

Animals undergoing treatment will most likely display a negative PCR result.

## ANALYTICAL SENSITIVITY

The PCRun® reaction can detect 10<sup>3</sup> copies of the target gene in pure DNA.

## REFERENCES

1. Comparison of PCR Assays for Diagnosis of Cutaneous Leishmaniasis. Esther Bensoussan,1 Abdelmajeed Nasereddin,1 Flory Jonas,2 Lionel F. Schnur,1 and Charles L. Jaffe1- J. Clin Microbiol, Apr. 2006, p. 1435-1439
2. Evaluation of PCR for diagnosis of visceral leishmaniasis. O F Osman, L Oskam, E E Zijlstra, N C Kroon, G J Schoone, E TKhalil, A M El-Hassan and P A Kage - J. Clin. Microbiol. 1997, 35(10):2454
3. Diagnosis of canine leishmaniasis: Conventional and molecular techniques using different tissues. Carla Maia, João Ramada, José M. Cristóvão, Luzia Gonçalves, Lenea Campino - The Veterinary Journal, Volume 179, Issue 1, January 2009, Pages 142-144



Manufacturer: Biogal Galed Labs. Acs. Ltd.  
TEL: +972 (0)4 9898605 FAX: +972 (0)4 9898690  
email: info@biogal.co.il www.biogal.co.il  
Kibbutz Galed, 1924000 - Israel