



**PCRRun<sup>®</sup>**  
by **Biogal** Galed Labs.

**PCRRun<sup>®</sup> brings  
gold-standard  
reliability to  
your clinic**



## Diagnostic Method



The PCRun kits, developed by Biogal Galed Labs, is a point-of-care molecular detection system for veterinary use, targeting infectious diseases in small animals. It uses an isothermal PCR method, which doesn't require a thermocycler, making it faster and more accessible for vet clinics.

## Benefits



### Point of Care

No extra equipment;  
works in clinics, hospitals,  
and labs



### Compact

Space-saving design for  
any lab or clinic



### Fast

Produces results in about  
an hour



### Accurate

High sensitivity,  
specificity, and scientific  
validation ensure reliability



### Affordable

Ideal for routine  
diagnostics without  
breaking the budget

















### User-Friendly

Minimal training required

## Applications

- ✓ **Symptomatic and Emergency Diagnostics** – Rapidly identify infectious pathogens in symptomatic animals or urgent clinical cases.
- ✓ **Outbreak Investigation & Containment** – Enable swift testing and response during disease outbreaks in kennels, shelters, and veterinary hospitals.
- ✓ **Pre-adoption or Breeding Screening** – Confirm infectious disease status before animal placement or reproduction.
- ✓ **Confirmatory Testing** – Support diagnosis when other tests yield inconclusive results or further validation is required.


# Analytical Sensitivity and Specificity of PCRun® Molecular Detection Kits

Product	Sample type	Sensitivity	Specificity	Study site
 <i>Babesia Species</i>	Whole blood	98.8%	97.5%	The Hebrew University of Jerusalem
 <i>Ehrlichia canis</i>	Whole blood	93.2%	100%	Department of Veterinary Medicine University of Bari
 <i>Anaplasma platys</i>	Whole blood	95%	100%	VBDDL North Carolina State University
 <i>Babesia gibsoni</i>	Whole blood	95%	100%	VBDDL North Carolina State University
 <i>Canine Distemper</i>	Nose swabs	100%	96.6%	CAVIDS Titer Testing Laboratory, University of Wisconsin
 <i>Pathogenic Leptospirosis</i>	Blood	92.6%	100%	Veterinary Diagnostic Laboratory, University of Wisconsin
 <i>Babesia canis</i>	Whole blood	91.2%	100%	Biogal Galed Labs
 <i>Leishmania infantum</i>	Whole blood/ lymph node	88%	100%	The Hebrew University of Jerusalem
 <i>Hepatozoon species</i> 	Whole blood	98.3%	100%	The Hebrew University of Jerusalem
 <i>Parvovirus</i> 	Plasma	100%	100%	CAVIDS Titer Testing Laboratory, University of Wisconsin
 <i>Feline Leukemia Virus</i>	Buffy coat/ plasma	93%	100%	Veterinary Diagnostic Laboratory, University of Wisconsin
 <i>Mycoplasma haemofelis</i>	Whole blood	95%	100%	Center for Companion Animal Studies, Colorado State University

All kits were validated using the gold standard Real-Time PCR method, with the exception of *Mycoplasma haemofelis*, which was evaluated using End-point PCR.


## STEP 1

## DNA EXTRACTION




## STEP 2

## DNA AMPLIFICATION



## STEP 3

## ANALYSIS OF RESULTS



Received: 26 March 2021 | Accepted: 20 January 2022

DOI: 10.1111/jvim.16373

### STANDARD ARTICLE

Journal of Veterinary Internal Medicine

Open Access



## Effect of sampling site on the diagnosis of canine parvovirus infection in dogs using polymerase chain reaction

Gilad Segev<sup>1</sup> | Tal Yaaran<sup>2</sup> | Sarah Maurice<sup>2</sup> | Gad Baneth<sup>1</sup>

Scan  
for full article



ISSN 1982-1263

<https://doi.org/10.22256/pubvet.v12n6a108.1-4>

## Detecção molecular de *Babesia canis vogeli* em cães da cidade de São Luís – MA, Brasil

Lygia Silva Galeno<sup>1\*</sup> , Brenda Fernanda Sodr  Moreno<sup>1</sup> , Andressa Mendes Alves<sup>1</sup> , Walkyria Concei  o Fonseca<sup>1</sup> , Clarissa Costa Dur es<sup>1</sup> , Douglas Marinho Abreu<sup>1</sup> , Italo Marcelo Reis da Silva<sup>1</sup> , Patr cia Thallyta Rocha Ferreira<sup>1</sup> , Daniel Praseres Chaves<sup>2</sup>

Scan  
for full article



Thomson et al. *Parasites & Vectors* (2018) 11:350  
<https://doi.org/10.1186/s13071-018-2914-5>

Parasites & Vectors

### RESEARCH

Open Access



## A new TaqMan method for the reliable diagnosis of *Ehrlichia* spp. in canine whole blood

Kirsty Thomson<sup>1\*</sup> , Tal Yaaran<sup>2</sup>, Alex Belshaw<sup>1</sup>, Lucia Curson<sup>1</sup>, Laurence Tisi<sup>1</sup>, Sarah Maurice<sup>2</sup> and Guy Kiddle<sup>1</sup>

Scan  
for full article



Biogal Galed Labs.

### HeatBlock

- ✓ **Compact and User-Friendly:** Perfect for in-clinic use.
- ✓ **Supports Two Critical PCRun® Steps:**
  - Step 1: DNA Extraction at 95°C for 5 minutes.
  - Step 2: Isothermal Amplification at 60°C for One Hour.
- ✓ **Dimensions:** 4 x 6 x 4.5 in. (10.5 x 15 x 11.5 cm)
- ✓ **Weight:** 2 lbs (1 kg)



### Qube4

- ✓ **Integrated Reader:** Designed for PCRun® amplification and analysis.
- ✓ **PC-Connected:** Enables user interaction, data analysis, and real-time results.
- ✓ **Multi-Sample Capability:** Processes four reactions simultaneously.
- ✓ **Compact Design:** Perfect for clinics and labs.
- ✓ **Dimensions:** 6.8 cm x 6.8 cm x 6.5 cm (2.68 in. x 2.68 in. x 2.56 in.)
- ✓ **Weight:** 180 g (0.4 lbs).



## CONTACT US



Biogal Galed Labs, Kibbutz Galed,  
1924000 Israel



+972 (0)4 9898605



[info@biogal.com](mailto:info@biogal.com)

For more information, visit our site  
[www.biogal.com](http://www.biogal.com)