VACCICheck[®]

by **Biogal** Galed Labs.

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JAVMA

A point-of-care dot blot ELISA assay for detection of protective antibody against canine adenovirus, canine parvovirus, and canine distemper virus is diagnostically accurate

Azizeh Egerer, BS; Zoie Schaefer, DVM; Laurie Larson, DVM*

Companion Animal Vaccines and ImmunoDiagnostic Service Laboratory, Department of Pathobiological Sciences, School of Veterinary Medicine, University of Wisconsin-Madison, Madison, WI

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Table 1—The sensitivity, specificity, and accuracy analysis of the index test for canine adenovirus (CAV), canine parvovirus (CPV), and canine distemper virus (CDV).

Virus	Sensitivity (95% CI)	Specificity (95% CI)	Accuracy (95% CI)
CAV n = 183 Pet dogs only	96.03% (91.55%–98.53%)	90.62% (74.98%–98.02%)	95.08% (90.87%–97.73%)
CPV n = 563 Pet and spf dogs	96.75% (94.39%–98.31%)	94.33% (90.08%–97.14%)	95.91% (93.93%–97.39%)
CDV n = 563 Pet and spf dogs	96.12% (93.69%–97.81%)	87.50% (81.69%–92.00%)	93.43% (91.05%–95.33%)

spf = Specific pathogen-free.

ANTIBODY TESTING VS VACCINATION
APPLICATIONS IN CLINICAL PRACTICE

RICHARD B. FORD.
PROFESSOR OF MEDICINE DIPLOMATE
ACVIM AND ACVPM (HON)

Escaneie para o artigo completo

Canine	
Adenovirus	In-clinic titer test results correlate well with gold standard testing (VN).
Distemper virus	In-clinic titer test results correlate well with gold standard testing (VN).
Parvovirus	In-clinic titer test results correlate well with gold standard testing (HI).

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Article

Prevalence of Serum Antibody Titers against Core Vaccine Antigens in Italian Dogs

Paola Dall'Ara ^{1,*}, Stefania Lauzi ¹, Jari Zambambieri ¹, Francesco Servida ², Laura Barbieri ³, Roberto Rosenthal ⁴, Lauretta Turin ¹, Elena Scarparo ¹ and Joel Filipe ¹

¹ Department of Veterinary Medicine and Animal Sciences (DIVAS), University of Milan, Via dell'Università 6, 26900 Lodi, LO, Italy

² Clinica Veterinaria Pegaso, Via Dante Alighieri 169, 22070 Rovello Porro, CO, Italy

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PAPER

Long-lived immunity to canine core vaccine antigens in UK dogs as assessed by an in-practice test kit

R. KILLEY*, C. MYNORS*, R. PEARCE*, A. NELL¹, A. PRENTIS² AND M. J. DAY^{3,1}

*Moreton Morrell College, Moreton Morrell, Warwickshire, CV35 9BL, UK

¹Medivet, Watford, Hertfordshire, WD24 7UY, UK

²Hyde Park Veterinary Centre, London, W2 2AE, UK

³School of Veterinary Sciences, University of Bristol, Langford, North Somerset, BS40 5DU, UK

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Article

Effect of Aging on the Immune Response to Core Vaccines in Senior and Geriatric Dogs

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Tutorial
no YouTube
CombCam

➤ Step 01

Use appropriate RFID
card with CombCam

➤ Step 02

Insert ImmunoComb
or VacciCheck Tooth

➤ Step 03

Read Result via
CombCam Cube or
PC Software

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