

ImmunoRun

Rapid Tests

Simple



Rapid



Accurate



Flexible



Biogal Gated Labs.

YOUR
DIAGNOSTIC
COMPANION

Lateral flow immunoassays play an important role in infectious disease diagnostics.

This technology allows veterinarians to provide an immediate and fast point-of-care diagnosis in the clinic or even in the comfort of the pet's home. The addition of our LFA new panels allows the differential diagnosis of infectious diseases that present similar clinical signs, leading to a rapid and more accurate treatment plan, increasing the success of a better prognosis.



Key advantages



Simple procedure with accurate results within 10 minutes.



High sensitivity, specificity, and accuracy.



Easy to use and interpret the results.













Economical test-kit, with convenient pack sizes of 5 or 50 tests.










Biogal Galed Labs.

YOUR
DIAGNOSTIC
COMPANION

Single Tests

Name	Sample	No. of Tests	Sensitivity	Specificity
 Immunodeficiency Virus Ab. Detection Kit	Whole blood, serum, or plasma	5 or 50	96.80%	98.90%
 Leukemia Virus Ag Detection Kit	Conjunctiva mucous membranes, nasal discharge, saliva, urine, serum or plasma	5 or 50	94.70%	99%
  Giardia Ag Detection Kit	Feces	5 or 50	100%	100%
 Parvovirus Ag Detection Kit	Whole blood, serum, or plasma	5 or 50	100%	100%
 Distemper Virus Ag Detection Kit	Feces	5 or 50	> 97.50%	> 97.50%
 Heartworm Ag Detection Kit (CDI)	Whole blood, serum, or plasma	5 or 50	94.40%	100%
 Leishmania Ab Detection Kit	Feces	5 or 50	95.60%	98%
  Corona Ag Detection Kit	Feces	5 or 50	99.99%	97.50%

Panels

Name	Sample	No. of Tests	Sensitivity	Specificity
  Crypto/Giardia Detection Kit	Feces	5 or 50	Crypto Ag 95%	Crypto Ag 97.50%
			Giardia Ag 91.89%	Giardia Ag 97.87%
 FeLV/FIV Detection Kit	Whole blood, serum, or plasma	5 or 50	FeLV Ag 94.59%	FeLV Ag 99.99%
			FIV Ab 91.67%	FIV Ab 97.83%
  Parvo/Corona Detection Kit	Feces	5 or 50	Parvo Ag 93.33%	Parvo Ag 99.99%
			Corona Ag 99.99%	Corona Ag 97.50%
  Parvo/Corona/Giardia Detection Kit	Feces	5 or 50	Parvo Ag 93.33%	Parvo Ag 99.99%
			Corona Ag 99.99%	Corona Ag 97.50%
			Giardia Ag 91.89%	Giardia Ag 97.87%



Scan
for more info



Complementary diagnostic solutions

VACCICheck[®]

Serology
Titer testing



Infectious Hepatitis
Parvovirus
Distemper Virus



Panleukopenia
Herpes Virus
Calici Virus



PCRun[®]

Molecular
Diagnostic PCR



Babesia Species
Pathogenic Leptospira

Ehrlichia canis
Anaplasma platys

Babesia canis
Babesia gibsoni

Parvovirus
Leishmania infantum
Distemper RNA



Parvovirus Panleukopaenia
Mycoplasma haemofelis
Leukemia Virus (DNA)

ImmunoComb

Serology
Dot Elisa



Parvo & Distemper Virus
Ehrlichia
Leptospirosis
Brucella



Feline *T. gondii* and *C. psittaci*
Feline Corona Virus (FIP)



Biogal Galed Labs.

YOUR
DIAGNOSTIC
COMPANION